



-Steam Distillation- -V1-

- Standard Operating Procedure -

Department - Process and Chemistry

Department Head - _____ **Date** _____

Provider - Gemstone Essential LLC

Purpose – To achieve a high quality fragrant essential oil of biomass sample for the purposes of compound qualification, product flavor/fragrance additive, and bioactive integration. Essential Oils are collected by creating conditions to quickly remove the Lighter-Than-Air components from a botanical sample, and removing the oil from atop its co-distilled hydrosol collection.

Document Control - Document Revision & Approval History					
Revision	Pages	Prepared By	Reviewed By	Approved By	Date
#1		Gemstone Essential	Gemstone Essential	Gemstone Essential	1/12/2019

Definitions –

Decarboxylation - Removal of the Carboxylic acids from compounds in the oleoresin

Essential Oil - Lighter-than-air and fragrant fractions of extracted oleoresin

Extraction – Treatment of organic matter with an organic or non-organic solvent for the purpose of concentrating valuable Drug and Non-Drug type molecules

Extraction Vessel- (EV) Container in which extraction protocols are conducted. An appropriate vessel that can handle required solvent, temperature, and pressure variables.

Extract Liquor – Post- extraction oil containing ethanol

Filtration – The use of filters, sieves, or separation to remove solid and semi-solid contaminants from a liquid collection

Keck clip – Plastic/ stainless steel clamp that fits snugly to specific scientific glass joint sizes to ensure separation does not occur

Oleoresin - Extracted Oil-bearing products from the Cannabis plant lineage

Terpenes – Individual Fragrant molecules that make up an Essential Oil

Volatiles - Fraction of low molecular weight compounds that have a BP below 120°C under vacuum.

Reaction - A process in which chemicals undergo a change of molecular structure

e.g. Decarboxylation of THC-Acid into THC

Equipment –

Boiling/ Liquid Trap Vessel (round bottom flask, jacketed vessel)

Heating equipment (heat mantle/ Steam generator, recirculating heater)

Material Vessel (Glass/ Stainless Steel biomass vessel)

Liquid Condenser (Glass/ Stainless)

Cold Heat Transfer Fluid Recirculator

Oil/water Separatory Funnel

Vapor Transfer Arm

Pipette


Essential Oil Storage Vessel (25-250mL)

Hydrosol Storage Vessel (1-5L)

Liquid transfer tubing

Procedure –

1. Break up Biomass until a 1-2 cm particle size is achieved throughout
2. Place a stainless steel mesh screen (10-25 lines per inch) at the bottom joint of the Material Vessel to ensure biomass is retained
3. Load Material Vessel with desired biomass. Use only moderate packing methods: (too tight of packing will cause steam channeling and cause a detriment to overall quality and yield)
4. Assemble Steam Distillation set up as shown below



- a. Fill Steam Generator with water
 - b. Attach water-out line of Cold Liquid Recirculator onto bottom hose barb of condenser
 - c. Attach water-in line of Cold Liquid Recirculator onto top hose barb of condenser
 - d. Attach Steam Generator steam hosing to Port 1 of Liquid Trap Vessel
 - e. Place Liquid Trap in Heat Mantle/ temp. appropriate vessel stand
 - f. Place Material Vessel on primary ground glass joint of Liquid Trap Vessel, secure with lab stand clamp and Keck Clip
 - g. Place large joint of Vapor Transfer Arm on top of Material Vessel joint
 - h. Place Condenser under small joint of Vapor Transfer Arm, Secure with Keck clip and lab stand clamp
 - i. Fill Oil/water Separatory funnel with Distilled water
 - j. Place Oil/Water Separatory funnel under bottom joint of Condenser. Secure with Keck clip and Lab stand clamp. Ensure stopcock valve is closed
 - k. Attach hose from hydrosol port of Oil/ Water Separatory Funnel and place in Hydrosol Collection Vessel
5. Set Cold Liquid Recirculator to 2-5°C and begin recirculation. If no Cold Liquid

- Recirculator is available, use a bucket of ice water and an aquarium pump
6. Begin Steam Distillation by initiating Steam Generator function
 - a. If no Steam Generator is available, Fill Liquid Trap under Material Vessel with distilled water, and set to 110°. Ensure flask is filled periodically.
 7. Allow steam to travel from Steam Generator through setup and into condenser
 8. If vapor is seen to be exiting Oil/Water Separatory Funnel, or if Oil/ Water Separatory Funnel is warm or hot to the touch, stop operation and check Condenser/ Cold Liquid Recirculator is operating properly.
 9. If liquid is trapped at the bottom of the Material Vessel, and appears to be soaking/ refluxing over material, momentarily pause Steam Generator operation and allow liquid to drain.
 10. Collect Essential Oil by removing Oil/ Water Separatory Funnel and using a pipette to transfer Oil into suitable Essential Oil Collection Vessel; Draining via Stopcock will cause loss of yield due to transfer loss on glass walls.
 11. Place air-tight cap on Essential Oil Collection Vessel and place in a freezer to freeze any remaining water/hydrosol. Filter through small glass funnel with clean a cotton ball into a new, clean Essential Oil Collection Vessel
 12. FOR BEST RESULTS:
 - a. Load only Distilled water into Steam Generator
 - b. After 20 minutes of operation: Collect multiple fractions of essential oil.
(Earlier fractions will be of higher quality)
 - i. Momentarily stop Steam Generator
 - ii. Remove Oil/ Water Separatory Funnel
 - iii. Use a pipette to drain off oil level which collects on top of the hydrosol collected within the Oil/ Water Separatory funnel
 - iv. Re-connect Oil/Water Separatory Funnel, ensure secure connection with Lab stand Clamp and Keck Clip
 - v. Resume operation
 - c. Use fresh, high quality biomass